

S100A10 Rabbit Pab - P60903

Item Number: CM0018891



Introduction

S100A10 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat, monkey S100A10. Ideal for phosphorylation signaling and calcium-binding protein research.

[Learn More](#)

Parameter	Value
Product Name	S100A10 Rabbit pAb
Remarks/Alias	42C, P11, p10, GP11, ANX2L, CAL1L, CLP11, Ca[1], ANX2LG, S100A10
Species	Human
Gene ID (Human)	6281
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-97 of human S100A10 (NP_002957.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	P60903
Protein Weight	11kDa
Shipping	Ice bag