


S100A5 Rabbit Pab - P33763

Item Number: CM0014105



Introduction

Rabbit polyclonal antibody against human S100A5 (S100D) validated for WB, IHC-P, and ELISA. Specific to human S100 calcium-binding protein A5, a metal ion (calcium, zinc, copper) binder. Ideal for research in signaling and metal homeostasis.

[Learn More](#)

Parameter	Value
Product Name	S100A5 Rabbit pAb
Remarks/Alias	S100D; S100A5
Species	Human
Gene ID (Human)	6276
Gene ID	6276
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-92 of human S100A5 (NP_002953.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P33763
Protein Weight	11kDa
Shipping	Ice bag