


S100B Mouse Mab For Wb - P04271

Item Number: CM0009316



Introduction

S100B Mouse Monoclonal Antibody validated for WB, IF-P, IHC-P, ELISA. Recognizes human, mouse, rat S100B (11 kDa), a calcium-binding protein highly expressed in astrocytes. Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	S100B Mouse mAb
Remarks/Alias	NEF; S100; S100-B; S100beta
Species	Human
GeneID (Human)	6285
GeneID	6285
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-92 of human S100B (NP_006263.1)
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04271
Protein Weight	11 kDa
Shipping	Ice bag