


## S100B Rabbit Mab For Wb, If, Ihc-P, Elisa - P04271

Item Number: CM0003649



### Introduction

S100B Rabbit Monoclonal Antibody, validated in WB, IF, IHC-P, ELISA. Detects human, mouse, rat S100B. Ideal for neuroscience, neurobiology, and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | S100B Rabbit mAb                                                                                                 |
| Remarks/Alias    | NEF; S100; S100-B; S100beta; S100B                                                                               |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 6285                                                                                                             |
| GeneID           | 6285                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-92 of human S100B (NP_006263.1). |
| Source           | Rabbit                                                                                                           |
| Category         | Monoclonal Antibodies                                                                                            |
| Application      | WB, IF-F, IF-P, IHC-P, ELISA                                                                                     |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | P04271                                                                                                           |
| Protein Weight   | 11kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |