

Saa3 Rabbit Pab - Saa3

Item Number: CM0014432

Introduction

SAA3 Rabbit pAb is a polyclonal antibody from rabbit, validated for WB and ELISA. It targets mouse SAA3 and cross-reacts with human, mouse, and rat SAA3, making it suitable for inflammation and acute-phase response research.

[Learn More](#)

Parameter	Value
Product Name	SAA3 Rabbit pAb
Remarks/Alias	I7R3; Saa-3; SAA3
Species	Mouse
Gene ID (Human)	6290
Gene ID (Mouse)	20210
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-123 of mouse SAA3P (NP_035445.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
Protein Weight	13kDa
Shipping	Ice bag