

Saa4 Rabbit Polyclonal Antibody For Wb, Elisa - P35542

Item Number: CM0020120

Introduction

SAA4 Rabbit Polyclonal Antibody (P35542) validated for WB and ELISA. This constitutively expressed serum amyloid A-4 protein antibody cross-reacts with human, mouse, and rat samples, ideal for acute phase research.

[Learn More](#)

Parameter	Value
Product Name	SAA4 Rabbit pAb
Remarks/Alias	CSAA; C-SAA; SAA4
Species	Human
Gene ID (Human)	6291
Gene ID	6291
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50 to the C-terminus of human SAA4 (NP_006503.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35542
Protein Weight	15kDa
Shipping	Ice bag