

Safb Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15424

Item Number: CM0030376

Introduction

SAFB rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA with human, mouse, rat cross-reactivity. Targets scaffold attachment factor B1 involved in transcription regulation, chromatin organization, and cell proliferation control.

[Learn More](#)

Parameter	Value
Product Name	SAFB Rabbit pAb
Remarks/Alias	HAP; HET; SAF-B; SAFB1; SAB-B1; SAF-B1; SAFB
Species	Human
GeneID (Human)	6294
GeneID	6294
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human SAFB (NP_001188267.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15424
Protein Weight	103kDa
Shipping	Ice bag