

Samd3 Rabbit Pab - Q8N6K7

Item Number: CM0020941

Introduction

Rabbit polyclonal antibody targeting SAMD3 (Sterile alpha motif domain-containing protein 3). Validated for WB and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for proteomics research and alternative splicing studies.

[Learn More](#)

Parameter	Value
Product Name	SAMD3 Rabbit pAb
Remarks/Alias	SAMD3
Species	Human
GeneID (Human)	154075
GeneID	154075
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human SAMD3 (NP_001017373.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N6K7
Protein Weight	61kDa
Shipping	Ice bag