

Sap18 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - 000422

Item Number: CM0025629



Introduction

Rabbit polyclonal antibody against human SAP18, validated for WB, IHC-P, and ELISA. Detects SAP18 in human, mouse, and rat samples. Ideal for research on transcriptional repression, splicing, and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	SAP18 Rabbit pAb
Remarks/Alias	SAP18P; 2HOR0202; SAP18
Species	Human
GeneID (human)	10284
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-172 of human SAP18 (NP_005861.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00422
Protein Weight	18kDa
Shipping	Ice bag