

Sars-Cov-2 Neutralizing Antibody - Sars-Cov-2

Item Number: CM0002104

Introduction

SARS-CoV-2 neutralizing monoclonal antibody, derived from recombinant protein immunogen, with endotoxin level <0.1 EU/μg by LAL method. Suitable for research applications involving SARS-CoV-2 neutralization studies. Shipped on dry ice.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 neutralizing antibody
Species	SARS-CoV-2
Immunogen	Recombinant protein
Source	SARS-CoV-2
Category	Monoclonal Antibodies
Application	< 0.1 EU/μg of the protein by LAL method.
Shipping	Dry ice