

Sars-Cov-2 Nsp1 Rabbit Pab - P0Dtc2

Item Number: CM0020073



Introduction

Rabbit polyclonal antibody targeting SARS-CoV-2 NSP1, validated for WB, IP, and ELISA. Reacts specifically with SARS-CoV-2; protein weight 141kDa. Ideal for viral infection research and assay development.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 NSP1 Rabbit pAb
GeneID	43740578
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of coronavirus NSP1 (YP_009742608.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTC2
Protein Weight	141kDa
Shipping	Ice bag