


Sars-Cov-2 Nsp6 Rabbit Pab For Wb And Elisa - P0Dtd1

Item Number: CM0024247



Introduction

Rabbit polyclonal antibody targeting SARS-CoV-2 non-structural protein 6 (NSP6). Validated for WB and ELISA, cross-reacts with SARS-CoV-2. Suitable for studying viral replication organelle formation and host-virus interactions.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 NSP6 Rabbit pAb
Gene ID	43740578
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 236-290 of coronavirus NSP6 (YP_009725302.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTD1
Protein Weight	794kDa
Shipping	Ice bag