

Sars-Cov-2 Nsp8 Rabbit Polyclonal Antibody For Wb, Elisa - P0Dtc2

Item Number: CM0020539



Introduction

SARS-CoV-2 NSP8 rabbit polyclonal antibody for WB and ELISA applications. Detects SARS-CoV-2 NSP8 protein; recombinant immunogen. Ideal for coronavirus research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 NSP8 Rabbit pAb
GeneID	43740578
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-198 of coronavirus NSP8 (YP_009742615.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTC2
Protein Weight	141kDa
Shipping	Ice bag