

Sars-Cov-2 Nucleoprotein Rabbit Pab - P0Dtc9

Item Number: CM0021822

Introduction

SARS-CoV-2 Nucleoprotein Rabbit Polyclonal Antibody (pAb) for WB, DB, IF/ICC, IP, ELISA. Raised in rabbit against a recombinant fusion protein covering amino acids 1-419 of coronavirus Nucleoprotein (P0DTC9). Specific to SARS-CoV-2.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 Nucleoprotein Rabbit pAb
Gene ID	43740575
Immunogen	Recombinant protein
SWISS	P0DTC9
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, DB, IF/ICC, IP, ELISA
Cross Reactivity	SARS-CoV-2
Shipping	Ice bag