

Sars-Cov-2 Orf3A Rabbit Pab - P0Dtc3

Item Number: CM0019825



Introduction

SARS-CoV-2 ORF3A Rabbit Polyclonal Antibody validated for WB and ELISA. Specific to SARS-CoV-2 ORF3a protein (P0DTC3). Ideal for studying viral egress, viroporin function, and host interaction.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 ORF3A Rabbit pAb
GeneID	43740569
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of coronavirus ORF3A (YP_009724391.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTC3
Protein Weight	31kDa
Shipping	Ice bag