

Sars-Cov-2 Orf9B Rabbit Polyclonal Antibody For Wb, Elisa - P0Dtd2

Item Number: CM0020350



Introduction

SARS-CoV-2 ORF9b rabbit polyclonal antibody validated for Western blot and ELISA. Targets the viral accessory protein ORF9b, which inhibits host innate immunity via mitochondrial association. Suitable for COVID-19 research.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 ORF9b Rabbit pAb
Species	SARS-CoV-2
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-97 of coronavirus ORF9b (P0DTD2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTD2
Protein Weight	11kDa
Shipping	Ice bag