

Sars-Cov-2 Spike Rbd (N501Y) Rabbit Mab For Wb, Elisa - P0Dtc2

Item Number: CM0002707



Introduction

High-specificity SARS-CoV-2 Spike RBD (N501Y) rabbit monoclonal antibody validated for WB and ELISA. Suitable for studying viral entry mechanisms and variant detection in research applications.

[Learn More](#)

Parameter	Specification
Product Name	SARS-CoV-2 Spike RBD (N501Y) Rabbit mAb
Remarks/Alias	spike glycoprotein; SARS-CoV-2 Spike RBD (N501Y)
Species	SARS-CoV-2
Gene ID (Human)	43740568
Gene ID	43740568
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 450-550 of coronavirus SARS-CoV-2 Spike RBD (N501Y) (YP_009724390.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTC2
Protein Weight	141KD
Shipping	Ice bag