

Sars-Cov-2 Spike S1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P0Dtc2

Item Number: CM0002525



Introduction

Rabbit monoclonal antibody against SARS-CoV-2 Spike S1 glycoprotein. Validated for Western blot, immunofluorescence, and ELISA. Targets a synthetic peptide within the S1 subunit (aa 200-300). For SARS-CoV-2 research.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 Spike S1 Rabbit mAb
Remarks/Alias	spike glycoprotein; SARS-CoV-2 Spike S1
Gene ID	43740568
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of SARS-CoV-2 Spike S1 (YP_009724390.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTC2
Protein Weight	141kDa
Shipping	Ice bag