

Sash1 Rabbit Polyclonal Antibody For Wb And Elisa - 094885

Item Number: CM0013589



Introduction

Rabbit polyclonal antibody against human SASH1, validated for WB and ELISA. Detects SASH1 in human and mouse samples. Ideal for studies on NF- κ B signaling, cell migration, and skin pigmentation.

[Learn More](#)

Parameter	Value
Product Name	SASH1 Rabbit pAb
Remarks/Alias	DUH; DUH1; CAPOK; SH3D6A; dj323M4.1; SASH1
Species	Human
Gene ID (Human)	23328
Gene ID	23328
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SASH1 (NP_056093.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O94885
Protein Weight	137kDa
Shipping	Ice bag