


Sbds Rabbit Pab For Wb, Elisa - Q9Y3A5

Item Number: CM0028834



Introduction

Rabbit polyclonal antibody against human SBDS, validated for WB and ELISA. Recognizes human, mouse, rat. Target involved in ribosome biogenesis and translation. Ideal for research on ribosome assembly and Shwachman-Diamond syndrome.

[Learn More](#)

Parameter	Value
Product Name	SBDS Rabbit pAb
Remarks/Alias	SDS; SDO1; SWDS; CGI-97; SBDS
Species	Human
Gene ID (Human)	51119
Gene ID	51119
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y3A5
Protein Weight	29kDa
Shipping	Ice bag