

Scg3 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Wxd2

Item Number: CM0030299

Introduction

Rabbit polyclonal antibody against human SCG3 (Secretogranin-3), validated for WB and ELISA. Detects human, mouse, and rat SCG3. Suitable for secretory granule biogenesis and angiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	SCG3 Rabbit pAb
Remarks/Alias	SGIII; SCG3
Species	Human
Gene ID (Human)	29106
Gene ID	29106
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WXD2
Protein Weight	53kDa
Shipping	Ice bag