


Scn2B Rabbit Pab For Wb, Elisa - O60939

Item Number: CM0024758



Introduction

Rabbit polyclonal antibody against SCN2B (sodium channel regulatory subunit beta-2) validated for Western blot and ELISA. Ideal for studying voltage-gated sodium channel modulation in excitable cells. Human cross-reactivity.

[Learn More](#)

Parameter	Value
Product Name	SCN2B Rabbit pAb
Remarks/Alias	ATFB14; SCN2B
Species	Human
GeneID (Human)	6327
GeneID	6327
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-160 of human SCN2B (NP_004579.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O60939
Protein Weight	24kDa
Shipping	Ice bag