

Scn4B Rabbit Polyclonal Antibody - Scn4B

Item Number: CM0021082

Introduction

SCN4B Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human, mouse, and rat SCN4B protein (~25kDa). Recombinant immunogen. Useful for neuroscience and cardiac research.

[Learn More](#)

Parameter	Value
Product Name	SCN4B Rabbit pAb
Remarks/Alias	Gm1471, SCN4B
Species	Mouse
Gene ID (Human)	6330
Gene ID (Mouse)	399548
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-228 of mouse SCN4B (Q7M729)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
Protein Weight	25kDa
Shipping	Ice bag