

Scrn1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q12765

Item Number: CM0019520

Introduction

High-quality rabbit polyclonal antibody targeting Secernin-1 (SCRN1, SES1) for Western blot, IF/ICC, and ELISA in human, mouse, rat samples. Ideal for studying exocytosis and mast cell biology.

[Learn More](#)

Parameter	Value
Product Name	SCRN1 Rabbit pAb
Remarks/Alias	SES1; SCRN1
Species	Human
Gene ID (Human)	9805
Gene ID	9805
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 195-414 of human SCRN1 (NP_055581.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12765
Protein Weight	46kDa
Shipping	Ice bag