

Scrn2 Rabbit Pab For Wb And Elisa - Q96Fv2

Item Number: CM0012826

Introduction

High-quality SCR2 rabbit polyclonal antibody validated for WB and ELISA. Suitable for detecting human and mouse Secernin-2, a 47 kDa phosphoprotein involved in multiple cellular processes.

[Learn More](#)

Parameter	Value
Product Name	SCR2 Rabbit pAb
Remarks/Alias	Ses2; SCR2
Species	Human
Gene ID (Human)	90507
Gene ID	90507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 176-425 of human SCR2 (NP_612364.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96FV2
Protein Weight	47 kDa
Shipping	Ice bag