


## Selenok Rabbit Pab For Wb And Elisa - Q9Y6D0

Item Number: CM0021902



### Introduction

SELENOK Rabbit pAb detects human and mouse SELENOK, a selenoprotein involved in calcium signaling, ERAD, and oxidative stress protection. Validated for WB and ELISA, suitable for immune and cardiac research.

[Learn More](#)

| Parameter        | Value                                                                                                |
|------------------|------------------------------------------------------------------------------------------------------|
| Product Name     | SELENOK Rabbit pAb                                                                                   |
| Remarks/Alias    | SELK; HSPC030; HSPC297; SELENOK                                                                      |
| Species          | Human                                                                                                |
| GeneID (human)   | 58515                                                                                                |
| GeneID           | 58515                                                                                                |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 1-94 of human SELENOK (NP_067060.2) |
| Source           | Rabbit                                                                                               |
| Category         | Polyclonal Antibodies                                                                                |
| Application      | WB, ELISA                                                                                            |
| Cross Reactivity | Mouse                                                                                                |
| SWISS            | Q9Y6D0                                                                                               |
| Protein Weight   | 11kDa                                                                                                |
| Shipping         | Ice bag                                                                                              |