


Serotonin Transporter Rabbit Pab For Wb, Ihc-P, Elisa - P31645

Item Number: CM0017474



Introduction

Human Serotonin transporter (SERT/SLC6A4) rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets N-terminal sequence. Suitable for neuroscience and serotonin signaling research.

[Learn More](#)

Parameter	Value
Product Name	Serotonin transporter Rabbit pAb
Remarks/Alias	HTT; 5HTT; OCD1; SERT; 5-HTT; SERT1; hSERT; 5-HTTLPR; Serotonin transporter
Species	Human
GeneID (human)	6532
GeneID	6532
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-84 of human Serotonin transporter (NP_001036.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31645
Protein Weight	70kDa
Shipping	Ice bag