

Serpinc1 Rabbit Pab - P01008

Item Number: CM0018384



Introduction

Rabbit polyclonal antibody against SERPINC1 (Antithrombin-III). Validated for WB and ELISA, with cross-reactivity for human, mouse, and rat. Ideal for blood coagulation research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	SERPINC1 Rabbit pAb
Remarks/Alias	AT3; AT3D; ATIII; THPH7; ATIII-R2; ATIII-T1; ATIII-T2; SERPINC1
Species	Human
GeneID (Human)	462
GeneID	462
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-259 of human SERPINC1 (NP_000479.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01008
Protein Weight	58kDa
Shipping	Ice bag