


Sesn2 Rabbit Pab For Wb, Elisa - P58004

Item Number: CM0025119



Introduction

SESN2 Rabbit pAb is a polyclonal antibody for WB and ELISA, targeting human SESN2 (cross-reacts mouse, rat). SESN2 is a leucine sensor regulating mTORC1 and oxidative stress; ideal for signaling and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	SESN2 Rabbit pAb
Remarks/Alias	HI95, SES2, SEST2, SESN2
Species	Human
GeneID (Human)	83667
GeneID	83667
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 211-480 of human SESN2 (NP_113647.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P58004
Protein Weight	54kDa
Shipping	Ice bag