

Setd6 Rabbit Monoclonal Antibody For Wb And Elisa - Q8Tbk2

Item Number: CM0006759



Introduction

Recombinant rabbit monoclonal antibody targeting human SETD6, validated for WB and ELISA. Detects human, mouse, and rat SETD6. Suitable for NF-kappa-B signaling and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	SETD6 Rabbit mAb
Remarks/Alias	SETD6
Species	Human
GeneID (human)	79918
GeneID	79918
Immunogen	Recombinant fusion protein containing amino acids 310-473 of human SETD6 (NP_001153777.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TBK2
Protein Weight	53kDa
Shipping	Ice bag