


Setd7 Rabbit Monoclonal Antibody For Wb And Elisa - Q8Wts6

Item Number: CM0011800



Introduction

Rabbit monoclonal antibody targeting the histone-lysine N-methyltransferase SETD7. Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Epigenetic research tool for histone methylation and transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	SETD7 Rabbit mAb
Remarks/Alias	KMT7; SET7; SET9; SET7/9; SETD7
Species	Human
Gene ID (Human)	80854
Gene ID	80854
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SETD7 (Q8WTS6).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WTS6
Protein Weight	41kDa
Shipping	Ice bag