


Sf3B2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13435

Item Number: CM0025145



Introduction

High-quality rabbit polyclonal antibody against human SF3B2, a key component of the spliceosome. Suitable for Western blot, IHC-P, IF/ICC, and ELISA; recognizes human, mouse, and rat SF3B2. Ideal for splicing research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	SF3B2 Rabbit pAb
Remarks/Alias	CFM; HFM; Cus1; SF3b1; SAP145; SF3B145; SF3b150; SF3B2
Species	Human
GeneID (Human)	10992
GeneID	10992
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 616-895 of human SF3B2 (NP_006833.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13435
Protein Weight	100kDa
Shipping	Ice bag