

Sfrp1 Rabbit Monoclonal Antibody For Wb, Elisa - Q8N474

Item Number: CM0010346

Introduction

SFRP1 Rabbit Monoclonal Antibody validated for Western blotting and ELISA, cross-reacting with human, mouse, and rat. Ideal for Wnt signaling and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	SFRP1 Rabbit mAb
Remarks/Alias	FRP; FRP1; FrzA; FRP-1; SARP2; SFRP1
Species	Human
GeneID (Human)	6422
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SFRP1 (Q8N474)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N474
Protein Weight	35kDa
Shipping	Ice bag