


## Sh2B1 Rabbit Pab For Ihc-P, If/Icc, Elisa - Q9Nrf2

Item Number: CM0030587



### Introduction

High-quality rabbit polyclonal antibody against human SH2B1 (PSM), validated for IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for studying JAK/STAT, insulin, leptin, and growth factor signaling pathways.

[Learn More](#)

| Parameter        | Value                                                                                                       |
|------------------|-------------------------------------------------------------------------------------------------------------|
| Product Name     | SH2B1 Rabbit pAb                                                                                            |
| Remarks/Alias    | PSM; SH2B; SH2B1                                                                                            |
| Species          | Human                                                                                                       |
| GeneID (Human)   | 25970                                                                                                       |
| GeneID           | 25970                                                                                                       |
| Immunogen        | Recombinant protein containing a sequence corresponding to amino acids 360-520 of human SH2B1 (NP_056318.2) |
| Source           | Rabbit                                                                                                      |
| Category         | Polyclonal Antibodies                                                                                       |
| Application      | IHC-P, IF/ICC, ELISA                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                           |
| SWISS            | Q9NRF2                                                                                                      |
| Protein Weight   | 79kDa                                                                                                       |
| Shipping         | Ice bag                                                                                                     |