


Shmt1 Rabbit Polyclonal Antibody For Wb, Ihc, If, Elisa - P34896

Item Number: CM0014861



Introduction

SHMT1 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat SHMT1 (53 kDa), a cytosolic serine hydroxymethyltransferase critical for one-carbon metabolism and folate cycle.

[Learn More](#)

Parameter	Value
Product Name	SHMT1 Rabbit pAb
Remarks/Alias	SHMT; CSHMT; SHMT1
Species	Human
GeneID (human)	6470
GeneID	6470
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 285-444 of human SHMT1 (NP_683718.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P34896
Protein Weight	53kDa
Shipping	Ice bag