

Shmt2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P34897

Item Number: CM0020679



Introduction

High-quality SHMT2 rabbit polyclonal antibody for WB, IHC-P, IF/ICC, IP, ELISA. Detects human, mouse, rat SHMT2 (~56 kDa). Ideal for one-carbon metabolism and mitochondrial disease research.

[Learn More](#)

Parameter	Value
Product Name	SHMT2 Rabbit pAb
Remarks/Alias	GLYA; SHMT; mSHMT; NEDCASB; HEL-S-51e; SHMT2
Species	Human
GeneID (Human)	6472
GeneID	6472
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 265-504 of human SHMT2 (NP_005403.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P34897
Protein Weight	56kDa
Shipping	Ice bag