

Siah1 Rabbit Mab - Q8luq4

Item Number: CM0004163



Introduction

SIAH1 rabbit monoclonal antibody validated for WB, IF-P, IHC-P, and ELISA. Detects human, mouse, and rat SIAH1 (31kDa). Ideal for ubiquitination, apoptosis, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	SIAH1 Rabbit mAb
Remarks/Alias	BURHAS; SIAH1A; SIAH1
Species	Human
GeneID (Human)	6477
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SIAH1 (NP_003022.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IUQ4
Protein Weight	31kDa
Shipping	Ice bag