

Siglec-15 Rabbit Pab - Q6Zmc9

Item Number: CM0019803



Introduction

Siglec-15 Rabbit Polyclonal Antibody for WB, IF-P, IHC-P, ELISA applications. Recognizes human and rat Siglec-15 (CD33L3), a 36 kDa immune receptor expressed in macrophage and dendritic cells.

[Learn More](#)

Parameter	Specification
Product Name	Siglec-15 Rabbit pAb
Remarks/Alias	CD33L3; HsT1361; SIGLEC-15; Siglec-15
Species	Human
Gene ID	284266
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-263 of human Siglec-15 (NP_998767.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q6ZMC9
Protein Weight	36kDa
Shipping	Ice bag