


Siglec1 Rabbit Pab - Q62230

Item Number: CM0023987



Introduction

Rabbit polyclonal antibody targeting mouse Siglec1 (Sialoadhesin, CD169) for WB and ELISA applications. Cross-reacts with rat. Serves as a key tool for studying macrophage adhesion, pathogen uptake, and immune response regulation.

[Learn More](#)

Parameter	Value
Product Name	Siglec1 Rabbit pAb
Remarks/Alias	Sn; Cd169; Siglec-1; Siglec1
Species	Mouse
GeneID (Human)	6614
GeneID	20612
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q62230
Protein Weight	183kDa
Shipping	Ice bag