


Sigmar1 Rabbit Pab For Wb, Ihc-P, Elisa - Q99720

Item Number: CM0026405



Introduction

Rabbit polyclonal antibody against human SIGMAR1 (sigma-1 receptor). Detects endogenous levels and works in WB, IHC-P, ELISA. Cross-reactivity with mouse and rat. Ideal for neuroscience, cell signaling, and neurodegeneration studies.

[Learn More](#)

Parameter	Value
Product Name	SIGMAR1 Rabbit pAb
Remarks/Alias	SRBP; ALS16; DSMA2; OPRS1; SR-BP; SIG-1R; SR-BP1; sigma1R; hSigmaR1; SIGMAR1
Species	Human
GeneID (human)	10280
GeneID	10280
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 114-223 of human SIGMAR1 (NP_005857.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99720
Protein Weight	25kDa
Shipping	Ice bag