

Sirt2 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Ixx6

Item Number: CM0018634



Introduction

SIRT2 rabbit polyclonal antibody validated for WB and ELISA, detects human, mouse, rat SIRT2. Ideal for studying SIRT2 deacetylase functions in cell cycle, transcription, metabolism, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	SIRT2 Rabbit pAb
Remarks/Alias	SIR2; SIR2L; SIR2L2; SIRT2
Species	Human
Gene ID (Human)	22933
Gene ID	22933
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 300-350 of human SIRT2 (NP_036369.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8IXJ6
Protein Weight	43kDa
Shipping	Ice bag