


Sirt6 Rabbit Polyclonal Antibody For Wb And Elisa - Q8N6T7

Item Number: CM0028895



Introduction

Rabbit polyclonal antibody targeting SIRT6, a NAD-dependent deacylase involved in DNA repair, metabolism, and aging. Validated for Western blot and ELISA; cross-reacts with human and mouse SIRT6.

[Learn More](#)

Parameter	Value
Product Name	SIRT6 Rabbit pAb
Remarks/Alias	Sir2I6; mSIRT6; 2810449N18Rik; SIRT6
Species	Mouse
GeneID (Human)	51548
GeneID (Mouse)	50721
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 200-294 of mouse SIRT6 (NP_001156902.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8N6T7
Protein Weight	36 kDa / 39 kDa
Shipping	Ice bag