

Skap2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - 075563

Item Number: CM0025774



Introduction

SKAP2 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets Src kinase-associated phosphoprotein 2 involved in B-cell adhesion and src signaling.

[Learn More](#)

Parameter	Value
Product Name	SKAP2 Rabbit pAb
Remarks/Alias	PRAP; RA70; SAPS; SCAP2; SKAP55R; SKAP-HOM; SKAP2
Species	Human
GeneID (human)	8935
GeneID	8935
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human SKAP2 (NP_003921.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075563
Protein Weight	41kDa
Shipping	Ice bag