

Skp1 Rabbit Polyclonal Antibody For Wb And Elisa - P63208

Item Number: CM0020992



Introduction

SKP1 rabbit polyclonal antibody (pAb) validated for WB, ELISA. Detects human and mouse SKP1, a core SCF ubiquitin ligase component (19 kDa). Suitable for cell cycle, signal transduction, and protein degradation research.

[Learn More](#)

Parameter	Value
Product Name	SKP1 Rabbit pAb
Remarks/Alias	OCP2, p19A, EMC19, SKP1A, OCP-II, TCEB1L, SKP1
Species	Human
Gene ID (Human)	6500
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-163 of human SKP1 (NP_733779.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P63208
Protein Weight	19kDa
Shipping	Ice bag