


Slc22A3/Oct3 Rabbit Pab For Wb, Elisa - 075751

Item Number: CM0019450



Introduction

Rabbit polyclonal antibody targeting human SLC22A3/OCT3 (EMT) for WB and ELISA, cross-reacts with mouse. Ideal for exploring organic cation transport, neurotransmission, and polyamine transport.

[Learn More](#)

Parameter	Value
Product Name	SLC22A3/OCT3 Rabbit pAb
Remarks/Alias	EMT; EMT3; OCT3; SLC22A3/OCT3
Species	Human
Gene ID (Human)	6581
Gene ID	6581
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 289-349 of human SLC22A3/OCT3 (NP_068812.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O75751
Protein Weight	61kDa
Shipping	Ice bag