

Slc22A4 Rabbit Pab - Q9H015

Item Number: CM0013768

Introduction

Anti-SLC22A4 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human, mouse, rat SLC22A4/ergothioneine transporter (OCTN1). Ideal for transporter and oxidative stress research.

[Learn More](#)

Parameter	Value
Product Name	SLC22A4 Rabbit pAb
Remarks/Alias	OCTN1; DFNB60; SLC22A4
Species	Human
GeneID (Human)	6583
GeneID	6583
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-141 of human SLC22A4 (NP_003050.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H015
Protein Weight	62kDa
Shipping	Ice bag