

Slc25A13 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Ujs0

Item Number: CM0025002



Introduction

High-quality rabbit polyclonal antibody against human SLC25A13 (citrin), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse. Suitable for mitochondrial transporter studies and liver disease research.

[Learn More](#)

Parameter	Value
Product Name	SLC25A13 Rabbit pAb
Remarks/Alias	CTLN2; NICCD; CITRIN; ARALAR2; SLC25A13
Species	Human
GeneID (Human)	10165
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SLC25A13 (NP_001153682.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UJS0
Protein Weight	74kDa
Shipping	Ice bag