

Slc25A6 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P12236

Item Number: CM0029074

Introduction

Rabbit polyclonal antibody targeting human SLC25A6 (ADP/ATP translocase 3, ANT3), suitable for WB, IF/ICC, and ELISA applications. Cross-reacts with mouse and rat. Ideal for mitochondrial transport and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	SLC25A6 Rabbit pAb
Remarks/Alias	ANT; AAC3; ANT3; ANT 2; ANT 3; ANT3Y; SLC25A6
Species	Human
Gene ID (Human)	293
Gene ID	293
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human SLC25A6 (NP_001627.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12236
Protein Weight	33kDa
Shipping	Ice bag