

Slc7A1 Rabbit Pab For Wb, Ihc-P, Elisa - P30825

Item Number: CM0019516

Introduction

High-quality anti-SLC7A1 rabbit polyclonal antibody for WB, IHC-P, ELISA. Validated on human, mouse, and rat. Research use for cationic amino acid transporter studies.

[Learn More](#)

Parameter	Value
Product Name	SLC7A1 Rabbit pAb
Remarks/Alias	ERR; ATRC1; CAT-1; HCAT1; REC1L; SLC7A1
Species	Human
GeneID	6541
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 530-629 of human SLC7A1 (NP_003036.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30825
Protein Weight	68kDa
Shipping	Ice bag