


Slco1B1 Rabbit Pab - Q9Y6L6

Item Number: CM0024501



Introduction

SLCO1B1 Rabbit Polyclonal Antibody validated for WB and ELISA. Targets human organic anion transporter SLCO1B1 (OATP-C), key in hepatic drug uptake and clearance. Ideal for liver transporter studies and drug interaction research.

[Learn More](#)

Parameter	Value
Product Name	SLCO1B1 Rabbit pAb
Remarks/Alias	LST1; HBLRR; LST-1; OATP2; OATPC; OATP-C; OATP1B1; SLC21A6; SLCO1B1
Species	Human
Gene ID (Human)	10599
Gene ID	10599
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 597-691 of human SLCO1B1 (NP_006437.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y6L6
Protein Weight	76kDa
Shipping	Ice bag